

## SCHEDA TECNICA

### LIFTNES®

#### Reusable Patches Neck and Décolleté with case

**LIFTNES® Reusable Patches** soft, thin and ergonomic, they adhere firmly to the skin, ensuring comfort and ease of wear. Designed to be reusable for a sustainable, low-waste solution, LIFTNES® Reusable Patches come in a practical, waterproof case with an airtight closure to protect them from humidity. The case is made of E.V.A., a lightweight, soft-touch material widely used in biomedical applications for its durability and reliable performance.

**HOW TO USE** Apply the patches to the targeted area for approximately one hour. After use, wash with mild soap, rinse with lukewarm water and gently pat dry. The patches can also be stored in the refrigerator.

#### TECHNICAL CHARACTERISTICS:

**Material:** E.V.A. (Ethylene Vinyl Acetate) case

The SKINIUS® case is made of E.V.A., a thermoplastic elastomer with flexible, resilient, rubber-like properties and a soft-touch finish. This material is widely used in biomedical applications for its durability and resistance to external agents. Lightweight, waterproof and resistant to humidity, E.V.A. combines comfort and performance and is practical, durable and hygienic for everyday use. Its elastic memory helps the case maintain its shape over time, even with frequent handling.

**Size:** 28x20 cm

**Color:** material E.V.A. with white border

**Customization:** gray logo print

#### TECHNICAL CHARACTERISTICS:

**Material:** silicon

**Color:** white

**Customization:** gray logo print

